

Ultra High Performance Liquid Chromatograph

# Nexera series



Mobile Phase Monitoring

Auto-Diagnostics


FlowPilot

Maximizing Throughput

Unattended Operation

Integrated Lab

Compact Design



© 2015 Waters Corporation

# EXPERIENCE NEW BENCHMARKS

- A New Benchmark of Intelligence
- A New Benchmark of Efficiency
- A New Benchmark of Design

Shimadzu has long been advancing the analytical performance of HPLC systems. At the same time, we recognize that overall efficiency depends not only on the performance of one instrument, but on the management of all devices within a lab. This realization leads us to now, a time in which AI capabilities have been incorporated to allow devices to detect and resolve issues automatically. In addition, lab management has been integrated using the Internet of Things (IoT) and device networking, making it simple to review the status of instruments and optimize resource allocation.

Building upon 40 years of experience in LC technology, the Nexera series is a family of UHPLC systems that marries these AI and IoT enhancements to set new industry standards in terms of intelligence, efficiency, and design.

# Nexera™ series

## Key Features

### Analytical Intelligence ANALYTICAL INTELLIGENCE

- Automated support functions utilizing digital technology, such as M2M, IoT, and Artificial Intelligence (AI), that enable higher productivity and maximum reliability.
- Allows a system to monitor and diagnose itself, handle any issues during data acquisition without user input, and automatically behave as if it were operated by an expert.
- Supports the acquisition of high quality, reproducible data regardless of an operator's skill level for both routine and demanding applications.





## Mobile Phase Monitor MPM-40

- Reservoir tray weight sensors monitor the remaining mobile phase in real time. [P. 7](#)

## System Controller SCL-40, CBM-40

- Supports remote monitoring via a smart device.
- Mobile phases can be purged and baselines checked easily via the touch panel.

## UV-VIS Absorbance Detector SPD-40 / SPD-40V Photodiode Array Detector SPD-M40

- Temperature control improves performance stability. [P. 13](#)
- Analytical data is linked to information about consumables to ensure traceability.

## Solvent Delivery Unit LC-40 Series

- Auto-diagnostics to detect problems during analysis, and an auto-recovery function. [P. 6](#)
- Reduces space requirements with a dual solvent delivery system. [P. 12](#)

## Autosampler SIL-40 Series PLATE CHANGER

- The injection speed is twice as fast as previous models, shortening multi-analyte processing times. [P. 10](#)
- Can perform continuous analysis on up to 44 microtiter plates. [P. 11](#)
- High reproducibility and ultra-low carryover for micro-volume injections. [P. 13](#)
- Automated sample preparation functions such as diluting samples, adding internal standards, and performing derivatization reactions, reduce labor.

## Column Oven CTO-40 Series

- Slim-type column oven with half the width of the previous model. [P. 12](#)
- Easy column attach/detach mechanism prevents peak broadening. (Nexlock™)
- Automatic management of column information such as number of injections and pressure during analysis.
- Monitoring of column degradation through usage history (when used with LabSolutions™).

# A New Benchmark of Intelligence

## Maximizing Reliability, Minimizing Down Time



### Fully Unattended Operation from Startup to Shutdown


The Nexera can be set to start up at a specified time, so that it can complete auto-purge, equilibration and baseline checks in advance, and be ready for analysis as soon as you arrive at the lab. Moreover, the system can be set up in advance to run without user intervention all the way from startup through analysis to shutdown.

You can view the status and predicted analysis completion time for multiple systems from any location via a smart device. None of these features requires any special software.

### Auto-Diagnostics and Recovery


In rare cases, air bubbles can form in the mobile phase and cause problems if inhaled into the pump. The Nexera has the ability to monitor baseline changes and pressure fluctuations to check for abnormalities.

When it detects an unusual fluctuation, it can automatically pause the analysis, purge the flow path, and restart analysis once it has confirmed recovery to normal pressure.





#### Auto-Diagnostics for Trapped Bubbles (Patent Pending)

If an air bubble becomes trapped in the pump head, it causes an abrupt drop in pressure, after which periodic pressure changes (pulsations) will occur. This can be detected using Shimadzu's proprietary bubble detection algorithm which assesses these distinctive pressure fluctuations caused by trapped bubbles to distinguish them from any expected changes in pressure.



## FlowPilot Protects Columns

HPLC columns can be damaged by sudden pump starts and stops or extreme gradient changes. The Nexera automatically uses FlowPilot (Smart Flow Control) to increase the flow rate gradually to the set point. There is no need to create startup protocols for each analysis.



### FlowPilot (Patent Pending)

The pump controls the flow rate according to oven temperature. The flow rate is gradually increased up to half of the set value. (1) The flow rate is kept constant until the oven is ready. (2) The flow rate is gradually increased up to the set value. (3)



## Mobile Phase Levels Measured in Real Time

Reservoir tray weight sensors (optional) can be used to monitor the volume of mobile phase or autosampler rinse solution in up to twelve containers. The containers can also be checked remotely from a smart device.

You will no longer need to worry about running out of mobile phase mid-analysis, because the device will notify you before starting the run if the volume remaining is too low.

### Comparing Solvent Volume Measurements with Predicted Consumption (Patent Pending)

The required mobile phase volume is predicted at the start of analysis from the weight sensor's measurement data, the instrument configuration information (the mobile phase used and the rinsing liquid configuration), and the analysis method. If a shortage occurs, the user is notified.



# A New Benchmark of Intelligence

## Remote Monitoring and Integrated Lab Management




### Reliably Track Column Usage\* ANALYTICAL INTELLIGENCE

Columns used in analysis need to be traceable to investigate any influence they may have on results. With the Column Management Device (optional), each column is recognized by an individual ID, and stats such as pressure and injection number are automatically recorded in the database and linked to each measurement data file. The device can be linked to all columns regardless of their brand or type.

The condition of a column can be checked directly with the column management software, so that you can start analysis without having to test column performance in advance.

\*To make use of this feature, a CMD compatible with the Nexera series and LabSolutions DB/CS are required.



#### Column Management Device (CMD) (Patent Pending)

The CMD (optional) compares the column information with the requirements of the analysis method and only begins analysis if they match. This device prevents mistakes in the choice of columns while at the same time storing the usage history of each column.



#### Step 1

Manage all your columns from the database, which can display usage status, injection count, etc. on one screen. Information about the column used is automatically linked to the analytical results.

#### Step 2


Deterioration of columns can be checked visually by examining peak shapes. Each compound listed in the compound table can be directly overlaid on the chromatogram.



## Check Chromatograms Online


View the operating status of instruments from a web browser using a smart device. This allows you to confirm parameters such as oven temperature and pump status directly from the web, and to monitor chromatograms in real time without returning to the lab.

In addition, SHIMADZU LabTotal™ Smart Service Net (optional) saves operating data from your instruments in a cloud server. Error information, including the date and time the errors occurred, can be sent via email.



Check the operating status of your LC systems

Select one system




Monitor chromatograms currently being taken




## Allocate Resources Efficiently

Manage the overall operation of your lab with SHIMADZU LabTotal Smart Service Net (optional). Review and compare instrument usage to maximize available analysis time, and manage shared consumables through a common system. (Patent Pending)



View graphs of monthly operation status



Zoom in further to view the daily operation status of each instrument



# A New Benchmark of Efficiency

## Automating Workflow, Maximizing Throughput





### Analysis Cycle Time Less Than 10 Seconds

The SIL-40 autosampler can process the entire injection cycle time in as little as seven seconds, twice as fast as the previous model. In addition, continuous analysis can be carried out on up to 44 MTPs (using 3 PLATE CHANGERS). Together these features dramatically increase analysis throughput.




Pharmacokinetic analysis requires not only speed but also high reliability at low concentrations. With its ultra-fast injection and ultra-low carryover, the SIL-40 autosampler delivers high reproducibility and reliability, even during an ultra-fast 30-second analysis.




## Quick and Reliable Mobile Phase Blending

The Nexera can automatically blend mobile phases at any set ratio. This speeds up the preparation of buffer solutions and the dilution of solvents, and can be used to easily prepare the exact amounts required for analysis, reducing waste as well as labor.



## Simplified Method Development

The Nexera Method Scouting System is capable of automatically switching between combinations of up to 6 columns and 8 mobile phases. It can run unattended, maximizing available uptime. In addition, Method Scouting Solution Software can automatically create methods with different columns, mobile phases, and gradient conditions from a single base method.





## Remote Instrument Maintenance Ensures Stable Operation



Each Nexera component automatically transfers the status of consumable parts, the traceability information of consumable parts after shipment from the factory, and various logs and error information to the SHIMADZU LabTotal Smart Service Net cloud system. Using this information, the system recommends the timing for consumable part replacement. Problems can be diagnosed remotely by Shimadzu service engineers. This means the system can be maintained in optimal condition and equipment management costs can also be reduced.

View the usage information for consumable parts in every device at a glance



Receive notifications before the replacement date

- Easy to replace consumables at the appropriate time
- Improved lab efficiency



Lab manager

- Balance maintenance schedule remotely
- Easy to determine appropriate maintenance



Service engineer



Start of use    Current status    Replacement  
Usage display for consumable pump parts

# A New Benchmark of Design

## Compact and Inventive



### Space-Saving Design


The Nexera frees up bench space with a compact design two thirds the size of Shimadzu's previous model.

### Energy-Saving Standby Mode



The Nexera uses over 80% less electricity when in standby mode, significantly reducing running costs and supporting an environmentally-friendly lab.

### Dual Injection Enables Simultaneous Analysis

Injection ports for two separate flow paths can be installed, allowing two different types of analysis (such as analysis of amino acids, organic acids or vitamins) to be performed using one system.



reddot award 2019 winner




#### Dual Injection (Patent Pending)


A second injection port and injection valve have been added. As a result, the sample can be injected independently into two analysis flow lines with a single needle, enabling the consolidation of a two-system analysis within a single system.

## High-Sensitivity Impurity Analysis

The SPD-M40 detector achieves an extremely high level of sensitivity and linearity (up to 2.5 AU). This allows quantitation of very low concentration impurities even in high-concentration samples. The UV cut-off filter installed in the detector prevents sample degradation due to UV light, helping to maintain good linearity at low concentrations.



Example: Ketoprofen analysis reveals impurities of area% < 0.2%



Detector linearity for ketoprofen

### Reduced Stray Light and Improved Measurement Circuits (Patent Pending)



With the SPD-M40, the measurement circuits have been improved, and the amount of stray light in the optical system, a cause for linearity errors, has been reduced (to less than 1/3 of the usual amount). As a result, the absorbance linearity has been improved to at least 2.5 AU.

## Ultra-low Carryover

The Nexera boasts ultra-low carryover, even on a high-sensitivity LC/MS/MS. This reduces time spent on rinsing, resulting in a shorter overall analysis time.

### Elimination of Dead Volume (Patent Pending)

Dead volume in the flow line will lead to carryover from sample residue. In the SIL-40 series, the sample loop and needle are connected linearly, eliminating dead volume, and reducing carryover.




Rinse Conditions: Normal

After analysis of Chlorhexidine, the carryover to a blank solution is negligible.

## Stable Baseline

Baseline fluctuations can affect peak area calculation, reducing the accuracy of quantitative results. The SPD-M40 photodiode array detector's "Advanced TC-Optics" function adjusts the temperature of the flow cell, lamp, and optical system to lessen the impact of external temperature changes. Noise and drift have also been reduced by 40% compared to the previous model.



### Advanced TC-Optics (SPD-M40) (Patent Pending)


With the SPD-M40, in addition to conventional cell temperature regulation, the lamp and spectrometer, which are responsible for thermal effects on absorbance, are also thermally regulated, further reducing baseline fluctuations caused by the light source.

## Faster Analysis Without Sacrificing Precision

It is important to use the most appropriate column in your LC system to achieve the highest efficiency and most accurate chromatograms. The Shim-pack series consists of a range of columns designed with the increased capabilities of the Nexera in mind. The superior ruggedness of Shim-pack columns ensures a long lifetime even with demanding sample matrices.

### Shim-pack Velox™

The superficially porous particle (SPP) technology allows for optimum separation and analysis times. The range of column types available enables you to choose the most appropriate column for each application.



	SP-C18	C18	Biphenyl	PFPP	HILIC
<b>USP classification</b>	L1	L1	L11	L43	L3
<b>Stationary phase</b>	Sterically protected C18	C18	Biphenyl	Pentafluorophenyl propyl	None
<b>Particle size (µm)</b>	1.8, 2.7, 5	1.8, 2.7, 5	1.8, 2.7, 5	1.8, 2.7, 5	2.7

\*To maximize column lifetime, columns with 1.8 µm particle size are recommended for use with a pressure under 80 MPa.

### Shim-pack™ XR-ODS II & III

Although the Shim-pack XR-ODS II has a particle size of 2.2 µm, it can be used up to a pressure of 60 MPa, making it appropriate for a wide range of analyses. The particle size of the XR-ODS III column has been reduced to 1.6 µm. With a maximum pressure of 100 MPa, it achieves high separation even with a short column length.

These columns are suitable for shortening analysis times while taking full advantage of the high separation power of UHPLC.

	XR-ODS II	XR-ODS III
<b>USP classification</b>	L1	L1
<b>Stationary phase</b>	C18	C18
<b>Particle size (µm)</b>	2.2	1.6

### Shim-pack™ GIS/GIST/GISS

The Shim-pack GIS/GIST/GISS series lineup includes columns with particle sizes ranging from 5 µm, appropriate for HPLCs, to 2 µm, appropriate for UHPLCs. With a variety of substrates available, they are ideal for method development.

**GIS Series:** HPLC columns packed with a high-purity silica gel material as a base. With high particle uniformity to secure the mobile phase path, these columns are ideal for low-pressure analysis.



**GIST Series:** An increase in the inertness of the silica results in improved peak shapes and ruggedness. Can be used in pH 1–10 environments. Easy to use for a wide range of chemicals and environments.

**GISS Series:** Adding to the advantages of the GIST series, a faster elution time has been achieved to provide even sharper peaks.


### Extract the Smallest of Shoulder Peaks

The i-PeakFinder tool can automatically recognize peaks even in complex chromatograms with high noise, and accounts for baseline drift for higher integration accuracy. Adjustable parameters in i-PeakFinder, such as peak detection threshold, allow the user to detect smaller peaks.


In addition, the i-PDeA II function can quantify peaks that cannot be completely separated by a column. Separate peaks by simply specifying the time and wavelength in LabSolutions.



The smaller peak is normally impossible to extract



i-PDeA II chromatogram identifies the two separate peaks




Spectrum analysis with i-PDeA II



### Consistent Peak Splitting




Manually setting the baseline for fused peaks is time-consuming, and the peaks may be split differently depending on the user. The i-PeakFinder has adjustable parameters for different analyses, allowing the consistent application of the best baseline in each situation.



### Dynamic Range Extension Function for Accurate Calibration

i-DReC is a new analytical method that significantly extends the dynamic range. It enables the analysis of high-concentration compounds without diluting them, and ensures a correct calibration curve.



Nexera, the Analytical Intelligence logo, Nexlock, LabSolutions, SHIMADZU LabTotal, Shim-pack, and Shim-pack Velox are trademarks of Shimadzu Corporation.



Shimadzu Corporation  
[www.shimadzu.com/an/](http://www.shimadzu.com/an/)

**For Research Use Only. Not for use in diagnostic procedures.**

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

Company names, products/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation, its subsidiaries or its affiliates, whether or not they are used with trademark symbol "TM" or "®".

Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The contents of this publication are provided to you "as is" without warranty of any kind, and are subject to change without notice. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication.